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G Guarantee Conditions

The manufacturer provides the end user with a 12-month guarantee, beginning from the date of invoice, covering all operational faults, or defects in materials or of manufacture.

In case of justified complaints the manufacturer or the Authorised Repair Centre will proceed with the repair or substitution of the product free of charge. To be able to make use of the free repair or substitution, there is an indispensable condition that, shipped together with the device, there is a proof of purchase document for the same device, which must include clearly legible references to the product, its serial number and date of purchase. The guarantee expires when any damages and/or their consequences may be attributed to unsuitable procedures or modifications of the product performed by third parties not authorised by the manufacturer: that is, if non-original spare parts or components are used.

The recognition of any other claims from any origin whatsoever, in particular any requests for indemnification for damages or interest are excluded.

The manufacturer may not be held responsible for damages, or injuries and their respective consequences derived:

- From excessive wear

- From connection of the instrument not compliant with the CE regulations.

- From improper tampering or maintenance performed by unauthorised personnel

- From use of non-original accessories or spare parts

- From lack of observance of the instructions for use, for assembly or maintenance, or from improper use of the product

- From unusual chemical, electrical or electronic influxes

- From faulty connections (air, water, electricity)

The guarantee does not cover the conductors such as the flexible "fibre optics", nor any other elements made of synthetic materials.

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GB (translation from the original in Italian)

Classification

Class IIa Medical device according to directive 93/42/EEC Directive and subsequent additions and amendments. This device has been manufactured in compliance with the IEC 60601-1.

Intended Use

Medical device, or a medical device accessory for professional use by Dentists, assisted by authorised personnel. The device is intended for use on Dental Units complying with the 93/42/EEC Directive and subsequent additions and amendments and installed in locations with an electrical plant system that complies with the IEC 60364-7-10.

Symbols

	Symbol to call the attention to other information found in the equipment User Manual.
	Symbol corresponding to "TYPE B PART APPLIED" according to IEC 60601-1. Indicates a protection grade against direct and indirect contacts.
	Part that may be sterilised in an autoclave
	Mark of Compliance with the 93/42/EEC Directive "MEDICAL DEVICES" and subsequent additions and amendments.
	Symbol for disposal pursuant to the 2002/95/EC, 2002/96/EC, 2003/108/EC Directives.

Precautions

*The device is not suitable for use with blends of inflammable anaesthetics and oxygen or nitrous oxide.

*Use of the device must respect the indications found in the 97050225 Instruction Manual supplied in the attached CD-ROM.

*The instrument is provided in a non-sterile State. Before use, see to its suitable cleaning, disinfection and/or sterilisation (as applicable).

*Use of electrically powered devices may interfere with pacemakers. In case of doubt regarding the treatment of patients with implanted cardiac stimulators consult the manufacturer's instructions for use.

*During use the device dust and fragments from the patient's oral cavity or from the same device may be projected into the environment (organic and inorganic dust and metal dust, fragments of the device or its tips and potentially infected biological material).

*Protect the patient, when possible, using a dental dam.

*Instruct the patient to breath through the nose in cases where the dam is not applicable.

*Medical personnel must put on suitable personal protection equipment.

*In case of visible damage or of the emission of unusual noises and/or vibrations or if overheating is noticed, do not use the device and call authorised technical service assistance.

*Before each use, verify the proper locking of the burs. If the burr does not fit correctly do NOT use them.

*Only use burns with a standard shank diameter (1.590 to 1.600 - ISO 1797-1)

*Do not use worn or damaged burns. Handle the burns carefully wearing protective gloves.

*The burr lock button overheat if pressed while the turbine is operating.

*Be very careful not to inadvertantly press it inside the patient's mouth.

*Use of the device must be in accordance with the CE Mark of compliance with 93/42/EEC Directive and subsequent additions and amendments, made with biocompatible materials complying with ISO 10933-1.

*Use only original spare parts and accessories, or those authorised by the manufacturer.

*Make sure that the surgical area cooling flow is present during operation. In the absence of cooling, suspend use of the device and contact the authorized technical service provider.

*LED radiation: do not stare into beam or view directly with optical instruments.

Class 1M LED equipment (400<λ<600nm), (IEC 60825-1:1993 + A1:1997 + A2:2001) (Pmax<1mW).

Assembling the burns

Simply press the middle of the button on the head to enable the burns to be inserted and removed. Press the button locks the burr. Check carefully that the burr is fully inserted and locked, pulling it with two fingers. Do not press the burr release mechanism if the burr is not completely still in order to avoid damaging the burr locking mechanism.

Perfekte Montage

Use of the device must be carried out with respect for the indications found in the specific 97050225 Instruction Manual supplied in the attached CD-ROM. The turbine must be lubricated at least once a day and before each time it is sterilized in the autoclave with the special "DAILY OIL" spray using the appropriate adapters. Follow the instructions on the can, making sure to apply the spray to two short jets. After lubrication, replace the turbine on the cord and activate it. Turn the burr lock button clockwise to lock the burr in place in order to avoid mechanical damage to the instrument, then spray the outside of the device with gauze or cotton. Regular use of the lubricant spray ensures the best operating conditions for the instrument. At least once a week, lubricate the "O-Rings" located on the quick coupling, using S1 silicone lubricant. Wearing disposable gloves, apply a layer of grease on the fingertips and lubricate using the fingers. The "O-Rings" are used to seal the quick coupling. They are used to supply the spray. It is recommended that the circuit be flushed with water at the end of each use. The manufacturer requires an annual check or overhaul of the device by authorised service assistance personnel.

Cleaning, Disinfection And Sterilisation

Before disinfection, remove the burr from the turbine and the turbine from the cord. To maintain the level of hygiene, after each use and within a short time, clean, disinfect and sterilize the turbine. Clean the spray nozzle spray both by cleaning the wire and by blowing in air using the syringe.

For cleaning and/or disinfection of the outside use gauze or cotton soiled in either water (700-1060 mbar) or alcohol (220 KPa, 2.2 bar) for 5 minutes (250 cycles tested - EN ISO 7785-1). Before each cycle in the autoclave, it is essential to lubricate the turbine using the "DAILY OIL" lubricant spray, applying two short jets. Reassemble the turbine on the cord and activate it in order to expel the excess lubricant. Dry the outside of the device with gauze or cotton and remove the quick coupling, if present, prior to sterilization.

*Do not sterilise the quick coupling.

*Do not use ultrasound cleaners.

*Do not immerse the turbine in disinfectant or sterilizing liquids.

*Never sterilise the handpiece in a dry heat steriliser.

*Do not leave the turbine in the autoclave at the end of the cycle.

*Periodically, check the autoclave according to the requirements of the manufacturer.

Sterilization temperatures above the limit may damage the device.

Reparations and overhauls

Installation, overhaul, calibration and repair of the device must be performed exclusively by technical service assistance personnel authorised by the manufacturer.

Intended Use

The method of application of the device is intermittent. The operational cycle indicated respects the following times:

Type of Instrument	Operation (minutes)	Pause (minutes)
Silent Power 2	5	5
Silent Power 2L	5	5
Silent Power 4L / 4L-K	5	5
Silent Power Miniature	5	5
Silent Power Gold	5	5

Environmental conditions for use:

- Ambient Temperature 10°-40° C

- Relative Humidity 30 - 75 %

- Atmospheric Pressure 700 - 1060 hPa (700 - 1060 mbar)

Environmental conditions permitted for transport and warehousing

- Ambient Temperature -20° +70° C

- Relative Humidity 10 - 90 %

- Atmospheric Pressure 500 - 1060 hPa (500 - 1060 mbar)

Scraping

The device must be carried out with respect for the laws in force and according to the individual national laws. The materials used for the manufacture are not harmful for human beings or animals that come into contact with or that are exposed to them.

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