

# Dr's Cam 7 moz

Dr's Cam 7 MOZ,  
newest technology you've never experienced before.



Two in One!  
Intraoral Camera and Mouse

## Features

- Built-in mouse
- Easy install (Real Plug & Play)
- Viewer included
- Angled handpiece design
- Software compatibility
- Slim head
- 330° Rotation
- OLED screen
- Windows compatibility
- Black & White function

## Dr's Cam 7 MOZ



## Features of Dr's Cam 7 MOZ



**Built-in touch mouse**  
Useable as an intraoral camera and mouse.



**Easy install (Real Plug & Play)**  
Driver not necessary. Simply plug and play.



**Dr's Viewer2**  
Image viewer included.



**Angled handpiece design**  
Familiar and comfortable design.



**Software compatibility**  
Keyboard short-cut data for using with other company viewer software is included.



**Slim head**  
Easy to get approach oral cavity.



**330° Rotation**  
330° rotating head helps to capture images more easily.



**OLED screen**  
Bright and clear OLED screen indicates various device states.



**Windows compatibility**  
Useable with any Windows OS.



**Black & White function**  
Take analog X-ray image to PC.

## Specification

|                    |                                                                       |
|--------------------|-----------------------------------------------------------------------|
| Size / Weight      | Approx. 230 x 27 x 25 mm<br>Approx. 120g<br>Cable length – Approx. 3M |
| Power supply       | 5V 500mA (USB power)                                                  |
| Resolution         | 1.3 Mega pixel (1280 x 1024)                                          |
| Sensor / Lens type | CMOS / 4 layers pin-hole lens                                         |

\* Specifications may change without prior notice.

## Minimum System Requirement

|           |                                     |
|-----------|-------------------------------------|
| PC System | 2.0 GHz CPU                         |
| RAM       | 1 GB RAM                            |
| OS        | Windows XP Service Pack 2 or better |
| USB port  | USB 2.0                             |

## Dr's Viewer2

Image Viewer included.

