

Vector – instruments, materials and accessories



COMPRESSED AIR
SUCTION
IMAGING
DENTAL CARE
HYGIENE

Materials



Paro Steribox silver

Sterilisation box with perforated base for sterilising and storing the Vector Paro handpiece and tool kits

2031-330-00



Scaler Steribox blue

Sterilisation box with perforated base for sterilising and storing the Vector Scaler handpiece and tool kits

2032-330-00



Vector Fluid polish

Polishing fluid with hydroxylapatite particles for smoothing and polishing; contents: 200 ml

CWZ510C2350



Vector/RinsEndo disinfectant

Ready-to-use disinfectant solution for cleaning and disinfecting hose connections, DGHM-listed; packaging size: 2.5 l bottle

CDZ501C6101



Vector cleaner

Cleaner for the Vector hose system
packaging size: 2.5 l bottle

CCA531A6150

Therapeutic concepts and instrumentation

The Vector system offers comprehensive possibilities for the treatment of periodontitis and periimplantitis as well as for professional teeth cleaning.

Periodontology and initial treatment



Tool kit Paro

Instrument case for 5 Paro instruments

2031-450-00



Paro Probe straight

Periodontal probe-equivalent instrument for the treatment of subgingival and supragingival labial, buccal, and oral surfaces

2030-151-01E Package contents: 3 pieces



Paro Lancet

Lancet-shaped instrument for the treatment of subgingival labial, buccal, and oral surfaces in pockets of up to 5 mm

2030-151-02E Package contents: 3 pieces



Paro Probe curved

Naber's probe-equivalent instrument for the treatment of interradicular spaces

2030-151-03E Package contents: 3 pieces



Paro Curette

Angled instrument for the treatment of subgingival and supragingival approximal surfaces

2030-151-04E Package contents: 3 pieces



Paro Probe Plus

Periodontal probe-equivalent instrument for the universal treatment of approximal, oral, and vestibular surfaces. Access depth of up to 5 mm in approximal area

2031-400-06E Package contents: 3 pieces

Periodontology, preservation, and periimplantitis treatment



Tool kit Recall/Implant

Instrument case for 3 Recall instruments

2031-460-00



Recall Probe straight CFK

Probe-shaped, semi-flexible instrument (carbon fibre material) for the treatment of subgingival and supragingival labial, buccal, and oral surfaces

2030-153-02E Package contents: 3 pieces



Recall Curette CFK

Curved, semi-flexible instrument (carbon fibre material) for the treatment of approximal surfaces


2030-153-05E Package contents: 3 pieces



Supra Probe flexible

Probe-shaped, highly flexible instrument for the gentle cleaning and rinsing of implant and tooth surfaces

2030-152-01E Package contents: 3 pieces



Periimplant soft

Flame-shaped, partly flexible, universally applicable instrument (carbon fibre material) for the gentle cleaning of sensitive implant surfaces.

2031-474-01E Package contents: 3 pieces



Periimplant hard

Conical-shaped instrument (carbon fibre material) for the reliable and efficient removal of dental plaque on implant surfaces.

2031-473-01E Package contents: 3 pieces

Prophylaxis



P4 Scaler instrument

Universal instrument for the removal of supragingival plaque

P4: 2032-414-00



P1, P2, P3 Scaler instruments


For the removal of subgingival plaque in deep gum pockets

P1: 2032-411-00 (straight)

P2: 2032-412-00 (curved to the right)

P3: 2032-413-00 (curved to the left)

Convince your patients!
We can provide free support



Vector marketing package

- Patient brochures with display
- Patient poster
- Prophylaxis cards

Vector DVD

- An application film explains how the instruments are used
- More documents and information about the Vector system

DÜRR DENTAL AG
Höpfigheimer Strasse 17
74321 Bietigheim-Bissingen
Germany
www.duerr.de
info@duerr.de

