

Whicam Story2

Enjoy Wireless!

Wireless Intraoral Camera with Mouse



Features

- High speed digital wireless
- Built-in mouse
- Angled handpiece design
- Big capacity battery
- Slim head / 330° Rotation
- OLED screen
- Compact size receiver
- Dental Viewer Software included



Whicam Story 2



Features of Whicam Story 2

-  **High speed digital wireless (800 Mbps/sec.)**
High speed digital wireless provides more clear live image without noise.
-  **Built-in mouse**
Useable as an intraoral camera and mouse.
-  **Angled handpiece design**
Familiar and convenient design.
-  **Big capacity battery**
Long time use with big capacity battery
-  **Slim head / 330° Rotation**
Easy to get approach anywhere oral cavity with slim head.
-  **OLED screen**
Bright and clear OLED screen indicates various device stats.
-  **Easy capture**
Can take image by capture button on the device without other equipment.
-  **Anti-fog**
It doesn't get fog on the lens.
-  **Software compatibility**
Keyboard short-cut data for using with other company viewer software is included.
* Most of software which can capture by keyboard short-cut is usable.
-  **Windows support**
Support Windows XP (32bit), Windows VISTA, Windows 7 (32bit & 64bit)
*Not support Windows XP 64 bit

Wireless receiver

USB power operation
It doesn't need other electric power supply. The receiver works with USB electric power only.

Compact design
Compact receiver saves space.

Unlimited receiver channel
Digital wireless type support unlimited channel. User can use many receivers with one camera.
* Analog type has just few channels. And should set channel adjustment.

Extension receiver holder

Specification

Wireless type	Digital wireless (800 Mbps/sec.)
Lens	Dual layer lens
Sensor type	CMOS / 1.3 Mega pixel
Battery	Lithium-ion (Samsung / LG cell)

Minimum System Requirement

PC System	2.0 Ghz CPU, 1 GB RAM
OS	Windows XP Service Pack 2 or better
USB port	USB 2.0

(Window XP 64 bit is not supported.)

* Internet access would need for some kinds of PC environments.

Dr's Viewer2

Image Viewer included.

