

Safety Data Sheet

according to regulation (EC) 1907/2006 (REACH)



Trade name: **PrintoSept-ID**

Page 1/5

Revised at: 02.09.2009

Print date: 22.02.2010

1. Identification of the substance/preparation and of the company/undertaking

Trade name: **PrintoSept-ID**

Recommended use: Cleaning and disinfecting solution, ready for use, free of aldehydes and phenols, for all dental impressions and worn dental prostheses as well as dental works to be soaked in.

Effect of the preparation: Disinfecting and cleansing agent

Manufacturer: ALPRO MEDICAL GMBH
Mooswiesenstr. 9
D-78112 St. Georgen (Germany)
Telephone: +49 7725 9392-0
Telefax: +49 7725 9392-91
Email: alpro@alpro-medical.de
Internet: www.alpro-medical.com

Technical contact point: +49 7725 9392-0

Emergency phone: +49 7725 9392-0 (business hours)
+49 761 19240 Poisoning-information-central, Freiburg (24h emergency call)

2. Hazards identification

This product does not require labelling in terms of preparation guideline 1999/45/EC.

3. Composition/information on ingredients

Chemical characterization: Mixture of substances listed below with nonhazardous additions in aqueous solution.

<i>Hazardous ingredients</i>	<i>CAS-No.</i>	<i>Weight-%</i>	<i>Code letter</i>	<i>R-phrases*</i>
Ethanol	64-17-5	15-30%	F	11
Dodecyldipropylenetriamine	2372-82-9	0,2-1%	C, N	22, 35, 50
Didecylmethylpolyoxethylammonium-propionate	94667-33-1	0,2-1%	C, N	22, 34, 50

*Full text of R-phrases: see under section 16

4. First aid measures

General information: Note precautions for self-protection.

On inhalation: Move affected person into fresh air and keep still and warm. In case of irritation of the airways, seek medical advice.

On contact with skin: Wash off immediately with soap and water. In case of skin reactions, seek medical advice.

On contact with the eyes: Rinse the eyes with water with the eyelids open for a sufficient length of time, then seek ophthalmological advice immediately.

If swallowed: Rinse out mouth and then drink plenty of water. Do not induce vomiting. Seek medical advice.

Information for the doctor: First aid, decontamination, symptomatic treatment.

Safety Data Sheet

according to regulation (EC) 1907/2006 (REACH)



Trade name: **PrintoSept-ID**

Page 2/5

Revised at: 02.09.2009

Print date: 22.02.2010

5. Fire-fighting measures

Suitable extinguishing media:	Water, foam, dry fire-extinguishing media, carbon dioxide
Extinguishing media which must not be used for safety reasons:	None known
Special exposure hazards arising from the substance or preparation itself, combustion products, resulting gases:	In case of heating or in case of fire toxic gases may be liberated.
Special protective equipment for firefighting:	Wear self-contained breathing apparatus. Do not inhale explosion and combustion gases.
Additional information:	Cool endangered containers with water spray jet and if possible remove off fire area.

6. Accidental release measures

Personal precautions:	Use personal protective clothing. High risk of slipping due to leakage/spillage of product.
Environmental precautions:	Contain spillage with absorbing material, e.g. with sand.
Methods for cleaning up/taking up:	Treat recovered material as prescribed in: Disposal considerations.

7. Handling and storage

Hints for safe handling:	Avoid contact with skin and eyes.
Hints for protection against fire and explosion:	Product is not flammable.
Further information:	Not necessary
Requirements for storage rooms and vessels:	Keep tightly closed and store only in original container.
Hints on storage assembly:	Not necessary
Further information on storage conditions:	Not necessary
Storage class:	LGK 10

8. Exposure controls / Personal protection

Constituents with occupational exposure limits resp. biological occupational exposure limits requiring monitoring:

Maximum occupational exposure concentration; see TRGS 900 (German regulation)

Ethanol	CAS-No.: 64-17-5	Air limit value: 500 ml/m ³ ; 960 mg/m ³
Overflow factors for short-term values:		Peak limitation-category 4

General health and safety measures:

- Do not eat, drink, smoke while working.
- Wash hands before breaks and on finishing work.
- Preventive skin protection using barrier ointment.
- Keep away of foods and drinks.
- Avoid contact with eyes, skin and clothing.
- Avoid inhalation of aerosols and vapours.

Safety Data Sheet

according to regulation (EC) 1907/2006 (REACH)

Trade name: **PrintoSept-ID**

Page 3/5

Revised at: 02.09.2009

Print date: 22.02.2010

Personal protection equipment:

Respiratory protection:	Not necessary if handled cautiously.
Hand protection:	E.g. natural rubber- or PVC-gloves are recommended.
Eye protection:	Tightly sealed goggles are recommended.
Protective clothing:	Not necessary if handled cautiously.

9. Physical and chemical properties

Form:	liquid
Colour:	yellow, clear
Odour:	of lemon
pH-value (undiluted):	10,5-11,5 (20°C)
Boiling point/Boiling range:	not determined
Melting point/Melting range:	not determined
Flash point (DIN, TRbF):	36°C
Flammability (solid, gases):	The product is not self-igniting.
Explosive properties:	The product is not explosive.
Lower explosion limit:	none
Upper explosion limit:	none
Ignition temperature:	not necessary
Oxidizing properties:	none
Vapour pressure:	not necessary at ...°C
Density:	0,965-0,975 g/cm ³ at 20°C
Bulk density:	not determinable at ...°C
Solubility in water:	soluble
in ethanol:	soluble
in hexan:	not soluble
Partition coefficient:	not necessary (n-octanol/water)
Viscosity (sort):	not necessary at ...°C
Electrical conductivity (undiluted):	300-500 µS/cm (20°C)
Refractive index nD:	1,3494-1,3509 (20°C)
Optical rotation αD:	not necessary

10. Stability and reactivity

Conditions to avoid:	No hazardous decomposition products if stored and handled as prescribed.
Materials to avoid:	No hazardous reactions with other substances known.
Hazardous decomposition products:	No decomposition if used as prescribed.

Safety Data Sheet

according to regulation (EC) 1907/2006 (REACH)



Trade name: **PrintoSept-ID**

Revised at: 02.09.2009

Print date: 22.02.2010

Page 4/5

11. Toxicological information

Acute toxicity:	No acute toxicity known.
Primary irritant effect on the skin:	No primary irritant effect on the skin known.
Primary irritant effect on the eye:	No primary irritant effect on the eye known.
Mutagenic potential:	No mutagenic potential known.
Sensitization:	No sensitizing effects known.

The product was classified on the basis of the calculation procedure of the preparation directive 1999/45/EC.

Based on manufacturer data of the main components the acute toxicity, the skin irritation, the mucous membrane irritation and the mutagenic potential of the preparation were evaluated. However, some data are not complete regarding particular main components. Nevertheless according to the experience of the manufacturer there are no other hazards expected than those which are already mentioned on the label.

12. Ecological information

Ecotoxicity:	Inhibition of degradation activity in activated sludge is not to be anticipated during correct introduction of low concentrations.
Persistence and degradation:	Ingredients are biodegradable or can be precipitated harmlessly.
Other adverse effects:	No other adverse effects known.
Further information:	Do not allow product to reach ground water or waters and do not empty undiluted product into sewers.

13. Disposal considerations

Disposal of the product:

In the EC member states there are no standardized regulations for the disposal of chemicals. In Germany the law of the cycle of economy and waste (KrW/AbfG) prescribes the principle of utilization and therefore a distinction must be made between "waste for utilization" and "waste for removal". Furthermore special cases – particularly in case of delivery – are also regulated by the federal states. Please contact the corresponding office (authority or waste disposal company) who will inform you about the utilization or removal.

Waste code no.*: 16 10 02

Disposal of uncleaned as well as completely emptied packings:

Disposal according to official regulations. Contaminated packings are to be treated like the substance.

Waste code no.*: 15 01 02

Further information:

If there is no official regulation, non-contaminated packings can be treated like domestic waste or can be recycled.

*Waste code no.: Waste incurred is classified following the code of the European waste register according to kind of waste and branch.

Safety Data Sheet

according to regulation (EC) 1907/2006 (REACH)



Trade name: **PrintoSept-ID**
Revised at: 02.09.2009
Print date: 22.02.2010

Page 5/5

14. Transport information

No hazardous material as defined by the transport regulations (ADR/RID, IMDG-Code, ICAO-TI/IATA-DGR)

15. Regulatory information

Labelling

Code letter: -
Hazard designation: -
contains: -
R-phrases: -
S-phrases: -

Special labelling of certain preparations (preparation directive 1999/45/EC annex V):

„Safety data sheet available for professional user on request.“

National regulations (Germany)

Restrictions of occupation:	None
Störfallverordnung:	-
Classification according to VbF:	Not subject to the VbF
Technische Anleitung Luft:	-
Water hazard class:	WGK 1 (slightly water hazardous); self-classification

16. Other information

Full text of R-phrases appearing in section 3: (This is not the classification of the product PrintoSept-ID!)

R11 Highly flammable.
R22 Harmful if swallowed.
R34 Causes burns.
R35 Causes severe burns.
R50 Very toxic to aquatic organisms.

Further information:

This information is based on our present knowledge and describes the safety measures which are to be taken for this product. However, this shall not constitute a guarantee for any product features.

Marking of changes:

A backslash in the left hand margin indicates an amendment from the previous version.