



dentac[®]
Protective Series

dentac

MAKE YOUR
Own Style



Index

Welcome to	5
<i>Moda</i> MASK	6-9
Surgical Mask 3S	10-11
Surgical Mask CLA	12-13
Surgical Mask IIR	14-15
Provisor Mask	16-17
UltraFace Shield	18-19
FoldaFace Shield	20-21
GoggleFace Shield	22-23
Provisor Goggle	24-25
Body Safe 2040	26-27
BodySafe 2045	28-29
Surgical Gowns	30-31
Anticore	32-33
Detro Air Spray	28-29
Detro Derm Lotion	30-31
Superox Antisepsis	32-33



dentac[®]
feel the innovation

Welcome to

dentac[®]
feel the innovation

With a 25 years of dental experience, we are here to serve you with our products designed and manufactured in Turkey, demanding to offer the experience to your needs globally. Since the day of the foundation till now, we are enriching our DENTAC brand's portfolio that have been keeping patented for 13 years with recently manufactured products. Having extreme help from our key opinion leaders and research and development department, DENTAC is already getting a real player in dental market that cover your needs from best Endodontic products line to Restorative products as well as other subsidiary products.

The company has two main lines in its portfolio as Endodontic and Restorative lines. With the patented technologies our files are manufactured to help to dentists in best possible ways to shorten the time of the operations and bring ease to dentists in all endodontic cases with subsidiary products. Having Restorative products completely manufactured in Turkey, we offer best quality and ease to our dentists globally in their restorative treatments.

Supplying all of our products complying with EN ISO 13485:2016, we give importance with ultimate care to our patients and dentists health. Keeping our belief alive in our beautiful country, a source of our motivation, we will be investing in manufacturing and we will be enriching our product portfolio with elegant top quality products

Your DENTAC team..

MAKE YOUR *Own Style*



SPECIAL DESIGN FOR YOU

Reflect your own style with limitless text and logo options.

SPECIAL DESIGN FOR YOUR COMPANY

You are able to advertise your institution safely.
Reflect your brand & logo.



FEEL SPECIAL WITH YOUR NAME

Become recognizable behind your mask.

MASK HOLDER

Mask holder for easy to use all day long.



COLOR OPTIONS FOR REFLECTING YOUR STYLE

Complete your style with tens of color options.

MAKE YOUR *Own Style*



" You can also design your own
Tmask by yourself from web page of

www.moda-mask.com "

* Available for Türkiye. Please ask for other countries.





T-Mask® Surgical 3S

Disposable Face Mask

Easy to use and has a required filter performance. Full automated **untouched** ultrasonic production. Provide exceptional comfort and protection for staff. Equipped with integrated plastic brace allowing perfect anatomic shaping, made from soft, non-irritating odourless material gentle to the skin free off glass fibres, latex and other allergenic substances. Layers water-repellent Non-Woven (PP) Fabric for protection from dust and moisture, Available with Blue, Green, Yellow, Red, Pink, Gray, Black, Beige, Dark Beige, Soft Blue, Soft Green, Soft Red, Soft Pink, Soft Yellow and White colors.



- 3 Layers
- Full ultrasonic shape
- Waterproof
- Ultra Soft
- Painless and adjustable ear hook with the lenght of 18cm x 10cm
- Soft Layer for Sensitive skin
- Non Steril
- Bacterial Filtration Efficiency (BFE) > %95
- In accordance with **EN 14683** and **EN 149** by accredited Lab.
- **Type I** according to **CE** Regulations
- **Level 1** according to **FDA** Regulations

- Exceptional Comfort
- Total Protection
- Hypo-Allergenic
- Comfortable Non-woven cotton lining



Advantages



Colors

○ White
● Blue
● Green
● Pink

● Beige
● Black
● Red
● Brown

● Yellow
● WaterBlue
● WaterGreen
● Light Pink

● WaterBeige
● Smoked
● LightRed
● LightBrown

● Light Yellow



Ref.No	Product Name	Ship. Type
DC2031-1	T-Mask Surgical 3S	1pcs/Box
DC2031-10	T-Mask Surgical 3S	10pcs/Box
DC2031-50	T-Mask Surgical 3S	50pcs/Box



T-Mask® Surgical CLA

Disposable Face Mask

Easy to use and has a high filter performance with %98 filtration. Full automated untouched ultrasonic production. Provide exceptional comfort and protection for dental, medical experts. Equipped with integrated plastic brace allowing perfect anatomic shaping, made from soft, non-irritating odourless material gentle to the skin free off glass fibres, latex and other allergenic substances. Layers water-repellent Non-Woven (PP) Fabric and Meltblown for protection from dust and moisture, available with Blue, Green, Yellow, Red, Beige, Black, Pink and White colors.



- 3 Layers with high level filtration
- Full ultrasonic shape for
- External Layer; waterproof
- Ultra Soft
- Painless and adjustable ear loop
- Soft Layer for Sensitive skin
- Mask sizes with 18x10
- Non Steril
- Bacterial Filtration Efficiency (BFE) > %98
- In accordance with **EN 14683** and **EN 149** by accredited Lab.
- **Type II** according to **CE** Regulations
- **Level 2** according to **FDA** Regulations

- Exceptional Comfort
- Total Protection with 98% filtration
- Hypo-Allergenic
- Comfortable Non-woven cotton lining



Advantages



Colors

○ White
● Blue
● Green
● Pink

● Beige
● Black
● Red
● Brown

● Yellow



Ref.No	Product Name	Ship. Type
DC2022-10	T-Mask Surgical CLA	10pcs/Box
DC2022-50	T-Mask Surgical CLA	50pcs/Box



T-Mask® Surgical IIR

Disposable Face Mask

Easy to use and has a ultimate filter performance with %99 filtration. Full automated untouched ultra-sonic production. Provide exceptional comfort and protection for dental, medical experts. Equipped with integrated plastic brace allowing perfect anatomic shaping, made from soft, non-irritating odour-less material gentle to the skin free off glass fibres, latex and other allergenic substances. Layers water-repellent Non-Woven (PP) Fabric and Meltblown for protection from dust and moisture, available with Blue, Green, Yellow, Red, Beige, Black, Pink and White colors.



%99
Filtration

- 3 Layers with high level filtration
- Full ultrasonic shape for
- External Layer; waterproof
- Ultra Soft
- Painless and adjustable ear loop
- Soft Layer for Sensitive skin
- Mask sizes with 18x10
- Non Steril
- Bacterial Filtration Efficiency (BFE) > %99
- In accordance with **EN 14683** and **EN 149** by accredited Lab.
- **Type IIR** according to **CE** Regulations
- **Level 3** according to **FDA** Regulations

- Exceptional Comfort
- Total Protection with 99% filtration
- Hypo-Allergenic
- Comfortable Non-woven cotton lining



Advantages



Colors

○ White
● Blue
● Green
● Pink

● Beige
● Black
● Red
● Brown

● Yellow



Ref.No	Product Name	Ship. Type
DC2023-10	T-Mask Surgical IIR	10pcs/Box
DC2023-50	T-Mask Surgical IIR	50pcs/Box



Provisor Mask®

Disposable Face Mask

Provisor Mask® With Gripper

Provisor mask offers high protection for nose, eyes and mouth during treatments for single use. Having the 100-micron visor thickness, the Provisor mask does not give mask sensation over the face and gives wide visibility to the user. The mask come forward with offering both mask and visor together and brings ease by composing them together. Simple design offers easy, quick use and %99 filtered protection.



- Easy to use
- Easy eyes, mounth and nose protection
- Clear and bright provisor
- Single use
- Available in both gripper and band attachment
- Made with high quality material
- Mask has %99 filtration capability.
- Only 100-micron visor thickness
- Bacterial Filtration Efficiency (BFE) > %99
- In accordance with **EN 14683** and **EN 149** by accredited Lab.
- **Type IIR** according to **CE** Regulations
- **Level 3** according to **FDA** Regulations

Purpose of use: For the protection of the eye, nose and mouth from the drop particles coming from outside and the bacteria carried in it.

- Offers protection from drop particles
- Available in two different design
- 100-micron visor thickness for wider vision



Advantages



Colors

- White
- Blue
- Green
- Pink

- Beige
- Black
- Red
- Brown

- Yellow



Ref.No	Product Name	Ship. Type
DC2035	Provisor Mask	25pcs/Box
DC2036	Provisor Mask With Gripper	25pcs/Box
DC2025	Mask Holder	



UltraFace® Shield

Full Face Protection

UltraFace Shield, which is suitable for daily use, is quite light. It is an ideal shield for protecting the face, eyes and nose outside. It is suitable for use in all daily activities. The rubber on the product makes it possible to use the product in teens and adults. The product which is very light can be easily put on and take off, and it never creates weight on the face. It does not prevent social communication during the period of use and it does not mist.



- Full Face Protection
- Suitable for Daily use (including outdoor)
- Cristal Clear
- Made with high quality materials
- Featherweight – Only 25gr.
- Full face cover – Eyes, Nose and Mouth
- With the rubber band it is usable for adults and children
- Quick and easy to use
- Mist and scratch resistant provide unrestricted visibility
- Ventiled foam designed for air flow and comfort
- Suitable for using with glasses
- Variable Colors
- Reusable
- Proprietary design
- 28x20 cm size
- Does not prevent social communication

It is suitable to be cleaned with a solution such as alcohol and reused.

The purpose of use : Should be for the protection of the face from the droplets and bacteria carried inside them.

- Aesthetic Design
- Perfect View
- Full Face Protection
- Durable Build
- Independent Use



Advantages



Colors

- White
- Black



Ref.No	Product Name	Ship. Type
DC2020	UltraFace Shield	Single Pack
DC2020/2	UltraFace Shield	Double Pack
DC2020/10	UltraFace Shield	Ten Pcs/Box



FoldaFace[®] Shield

Full Face Protection

FoldaFace, which is suitable for daily use, is quite light. Shield can be adjusted at different angles up to 180 degrees. Its easy design ensures constant face protection. The acetate of the product does not mist. It is suitable for cleaning and disinfection. Its design allows for easy put on and take off. Its rubberized structure is designed for everyone to use. In case of need, the shield can be shifted and removed from the face without needing to be removed from the head. Since it does not create pain and pressure in the head, it is suitable to be worn on the face for a long time.

- Suitable for Daily use (including outdoor)
- Made with high quality material
- 180° Foldable Shield
- Reusable frame and quick-change disposable shield.
- Suitable for use with glasses
- Durable for scratches
- Comes with 1 frame and 5 shields or 3 frame and 15 shields
- Full face protection
- Mist and scratch resistant provide unrestricted visibility
- Quick and easy to use
- Available in 4 colors
- Proprietary design
- Does not prevent social communication

It is suitable to be cleaned with a solution such as alcohol and reused.

The purpose of use : Should be for the protection of the face from the droplets and bacteria carried inside them.

- 180° Foldable shield
- Perfect View
- Full Face Protection
- Durable Build
- Suitable for using all day

Colors

● Turquoise

● Red

○ White


● Black



Advantages

CE	ISO 9001	ISO 13485	EN 166: 2001	FDA
----	----------	-----------	--------------	-----

Ref.No	Product Name	Ship. Type
DC2010	FoldaFace Shield	Single Pack
DC2010/3	FoldaFace Shield	Triple Pack



GoggleFace® Shield

Full Face Protection

The easy and simple design of the product provides quick facial protection. The shield of the product does not mist. It is suitable for cleaning and disinfection. Its design allows for easy put on and take off. Since its design is in the style of glasses, there is no need for any uncomfortable additional connections that will surround the head. It can be quickly take-off from the head when needed. Since it is worn in the form of glasses, it does not create pain and pressure in head. It is suitable to be worn on the face for a long time. Because of its design, it is quite aesthetic. It can be used especially by specialists without hesitation throughout the day.

- Goggle design
- Elegant look
- Suitable for Daily use
- Full face protection
- Made with high quality material
- User friendly design
- Durable for mist and scratches
- Ultra-light 16gr
- Comes with filmed 5 shields and 1 frame (or 25 shield and 5 frame)
- Adaptable for glasses with extra apparatus
- Durable Fog and scratches

It is suitable to be cleaned with a solution such as alcohol and reused.

The purpose of use: Should be for the protection of the eye from the droplets and bacteria carried inside them.

Packing: 1 Frame, 5 Shields / 5 Frames, 25 Shields

- Perfect View
- Film Covered Shield
- Full Face Protection
- Durable Build
- Suitable for using all day



Advantages



Colors

- Blue
- Green
- Pink
- White

- Transparent
- Turquoise
- Red



Ref.No	Product Name	Ship. Type
DC2030	GoggleFace Shield	Single Pack
DC2030mix	GoggleFace Shield	5 pc/Box



Provisor Goggle

Eyes Protection

Provisor goggles especially protect your eyes from small impacts, vapors, dust, gas and liquids. It has a silicone interior design for a long-time use. It is suitable for re-use when it is cleaned and disinfected. Thanks to the wide spectacle structure, the angle of vision is quite wide. It has an adjustable design for every head size. Unisex design allows it to be used by anyone. The product available in lightweight design do not create weight on the head and face. The product material in the hypo allergenic composition minimizes the allergen risk. The nose pads are extra soft so they do not cause any discomfort in the nose during long use.

- Leak-proof design
- Protects against liquids, dust, gases and vapors
- Adjustable head band
- Doesn't mist
- Provides complete protection against small impacts
- Can be used in interior environment
- Silicone interior design closes all around the eyes and eliminates the discomfort
- Suitable for all face types
- Consists of high-quality material
- It has wide angle of vision
- Lightweight design does not create pressure on the face

Purpose of use: For the protection of the eye from the drop particles coming from outside and the bacteria carried in it.

- Waterproof, steam proof and dustproof
- Complete eye protection against minor impacts
- Can be disinfected
- Hypo-allergenic
- Provides wide angle of vision



Advantages



Colors

- Blue
- Green
- Pink
- Transparent

● Turquoise



Ref.No	Product Name	Ship. Type
DC2030/1	Provior Goggle	Single Pack
DC2030/5	Provior Goggle	5pc/Box



BodySafe® 2040

Disposable Protective Overalls

Dentac BodySafe 2040 disposable coverall protects against infected agents, especially those that can be transmitted through droplets. Spunbond PP Nonwoven Fabric is used for long term use. It is not suitable for re-use. It provides comfortable movement thanks to its elastic waist, arm and foot design. Unisex design allows it to be used by anyone. For full protection, it is recommended to cover the open hand, foot and face periphery of the body with suitable materials such as; **Ultraface Shield, FoldaFace Shield, GoggleFace Shield, Provisor Goggle, T-Mask, Profeet Shoe Cover, Profeet Boot Cover.**

- Droplet-proof design
- Protects against droplets, dust and infected agents
- Elastic waist, arms, face circumference and ankle design
- Provides comfortable movement
- Easy to wear
- Breathable Spunbond PP Nonwoven fabric (53 gr/m²)
- Provides protection against chemical electrostatic pollution
- Light proof
- Latex-free
- Adhesive front zipper tape
- Hooded design, including chin
- Made of high quality material
- Interior design with sweat absorbing feature that does not irritate the skin
- Does not deposit ruins on working area
- Lightweight design does not cause discomfort in long-term use.

Purpose of use: For the protection of the body from the drop particles coming from outside and the bacteria carried in it.

Packing: Available in individual bag in White color.
Dentac UltraFace, gloves and single use Profeet Boot Galosh are not included in the package.

Sizes: S, M, L, XL, XXL



Advantages

- Air, droplet and dustproof
- Sealed Tape design
- Comfortable use
- It is not allergenic since it does not contain latex.



Ref.No	Product Name	Ship. Type
DC2040	BodySafe 2040	Single Pack
DC2050	Profeet Boot Galosh	Single Pack
DC2055	Profeet Shoe Galosh	Single Pack



BodySafe® 2045

Disposable Protective Overalls

Dentac BodySafe 2045 disposable coverall protects against infected agents, especially those that can be transmitted through droplets. Spunbond PP Nonwoven Fabric is used for long term use. There are air, water and gas tight special elastic bands at the seam joints. It is not suitable for re-use. It provides comfortable movement thanks to its elastic waist, arm and foot design. Unisex design allows it to be used by anyone. For full protection, it is recommended to cover the open hand, foot and face periphery of the body with suitable materials such as; **Ultraface Shield, FoldaFace Shield, GoggleFace Shield, Provisor Goggle, T-Mask, Profeet Shoe Cover, Profeet Boot Cover.**

- Droplet-proof design
- Special elastic banded design at the seam joints
- Protects against droplets, dust and infected agents
- Elastic waist, arms, face circumference and ankle design
- Provides comfortable movement
- User Friendly
- Breathable Spunbond PP Nonwoven Fabric (53 gr / m²)
- Provides protection against chemical electrostatic pollution
- Lightproof
- Does not contain latex
- Adhesive front zipper tape
- Hooded design, including chin
- Made of high-quality material
- Interior design with sweat absorbing feature that does not irritate the skin
- Does not leave fiber spillage in the working environment
- Lightweight design does not cause discomfort in long-term use.
- Latex Free

Purpose of use: For the protection of the body from the drop particles coming from outside and the bacteria carried in it.

Packing: Available in individual bag in White color. Dentac FoldaFace, gloves and single use Profeet Boot Galosh are not included in the package.

Sizes: S, M, L, XL, XXL



Advantages

- Waterproof and dustproof
- Sealed Tape design
- Comfortable use
- It is not allergenic since it does not contain latex.
- Comfortable front zipper design with protection tape



Ref.No	Product Name	Ship. Type
DC2045	BodySafe 2045	Single Pack



SMS Disposable Surgical Gowns

Body Protection

Dentac SMS Disposable Surgical Gown with Ribbon is made of high-quality 40gr/m2 SMS fabric and offers great comfort while using the products. It is suitable for daily use and use in hospitals and it can be easily worn on a daily dress. Available in 2 different colors with light blue and white, so it appeals to all medical requests. Optionally, the product has with and without ribbon models. The products can be used by doctors and patients. The belt over the product helps to tighten it to be fit by everybody. The product comes in plastic packaging with product label showing the products' features.

- SMS Fabric (Polypropylene Spunbond Nonwoven Fabric + Meltblown Nonwoven Fabric + Polypropylene Spunbond Nonwoven Fabric)
- Non Sterile (Autoclavable at 134o)
- Anti-Alcohol, Anti-Blood, Anti-Oil, Anti-bacteria
- Waterproof, Acidproof, Alkaliproof
- Made of 40 g/m2 SMS
- Arms are available in both ribbon and rubbery models
- Can be worn easily on the dress
- Available in Light blue and White
- Breathable, Flexible
- Available in XL size 150x120 cm (Can also be provided in customized sizes)
- Unisex Design

Purpose of use: For the protection of the body from the drop particles coming from outside and the bacteria carried in it.

- Durable
- With or without ribbon design
- Easy to wear
- Breathable & Comfortable
- High quality isolation



Colors

- Light Blue
- White



Advantages



Ref.No	Product Name	Ship. Type
DC2070	Surgical Gown	Single Pack
DC2071	Surgical Gown With Ribbon	Single Pack



Anticore

Body Protection

Anticore products that cover the nose, face and mouth completely protect you up to waist level. While the non-woven fabric provides air permeability, it provides liquid sealing. It is suitable for daily use and use in hospitals. It can be easily worn on a daily dress. Available in 4 different colors, so it appeals to all tastes. Because of its design, it is quite aesthetic. It can be used without hesitation throughout the day. The products in different sizes are suitable for children and adults. The kangaroo pocket on the product helps you store your personal products on the product without touching different areas from the outside. The product comes in plastic packaging. The product folds out of the kangaroo pocket, saving space in your closet.

- Nonwoven Fabric
- The arms and face of it are rubbery.
- Stitched with a Nano needle.
- Ultra-Light
- 45 gr / m² Fabric (or Customized Fabric)
- Can be worn easily on the dress
- Closes from head to waist
- Available in Navy Blue, Red, Black and White
- Breathable, liquid proof
- Presented with kangaroo pocket
- Available in S-M-L-XL
- Unisex Aesthetic Design

The purpose of use: Should be for the protection of the body and certain parts of the face from the droplets and bacteria carried inside them.

Packing Way: S-M-L-XL sizes in single packages

- Ultra-Light
- Waterproof
- Breathable
- Easy to wear
- Available in 4 different Colors
- Unisex Design



Advantages



Ref.No	Product Name	Ship. Type
DCAT2010	AntiCore	Single Pack



Detro Air Spray

Area Disinfection

- It is proper for disinfection of areas between 1-1000 m³.
- Easy to start and adjust the program to be applied through the keypad on.
- Auto printing action through it's integrated printer. It can also print 100 pcs of operation retrospectively.
- It is easy to operate the machine through the short cut keys and gives the chance to the user to save records for up to 100 rooms through "List Recording" function.
- It gives the chance to set the timing for startup.
- It can be controlled remotely from to distance of 50-100 meters through the optional RF controller.
- The operation can be paused and resumed just in case.
- Adjustable screen brightness.
- Passwords can be adjusted for safety.
- In order to be performed easy and effective pulverization the head is positioned by angel of 45 degrees.
- Visible and audible warning systems both in the beginning and at the end of operations.
- English language is available.
- Easy to hold through plastic grips.
- Particle filter keeps the air Access system clean.

Detrosan Air

- Stabilized hydrogen peroxide based disinfectant solution for
- automted airborne room disinfection bactericidal and yesticidal spectrum.
- Aldehyde free
- Leaves no trace, wetness or grassy film
- %99 biodegradable, ecofriendly formula
- Non-allergenic
- Bactericidal, fungicidal, virucidal, sporicidal
- Wide compatibility with surfaces and materials, noncorrosive

Composition:

- Stabilized hydrogen peroxide,
- Silver salt complexes

Action Mechanism

- Hyrdogen peroxide (H2O2) is converted quickly into hydroxyl radicals (OH-) those react immediately with the essential biomolecules (proteins, liquids, nucleic acids, etc.) inactivating them.
- Silver cations (Ag+), binding particularly with proteins, denaturing enzymes and interfering with DNA replications.
- The combined action of two active principles leads to cell inactivation to destruction of viruses, bacterias and spores.

Efficiency	Standards	Bactericidal	EN 13727	Fungicidal
EN 13624	Virucidal	EN 14476	Tuberculocidal	
EN 14348	Sporicidal	EN 13704		

It is used in ground and surface disinfection by pulverization method in critical areas such as operating room, intensive care, burn unit, infected patient rooms, laboratory and all general areas with Detro San Air solution.

Purpose of use: To disinfect all areas up to 1000 m3 (clinic, waiting rooms, stores etc) including air and all materials available inside.

Packing: Detro Air Spray device, power supply, remote control

Advantage: Offers affective disinfections in all areas.

CE	ISO 9001	ISO 22716		
EN 13624	EN 14348	EN 14476	EN 13704	EN 13727

Ref.No	Product Name	Ship. Type
DC2061	Detro Air Spray	
DC2061-S1	Detro San Air	1lt
DC2061-S20	Detro San Air	20lt



Detro Derm Lotion

Hand Antiseptic

- Alcohol based hand sanitizer
- Hand and skin sanitizer with alcohol affective in short time applied by spreading it to hand and skin.
- As there is no active materials cumulating over the skin in its formula, it can be easily tolerated by skin.
- Tested dermatologically.
-

Certificates : CE, ISO 9001, 22716

Advantage:

- Easy use due to distinctive dispenser.
- Ready to use

Composition: 2-propanol (%70V/V), glycerin, excipients

Microbiological Affect: Bactericidal (Including Tbc), fungicidal, virucidal, tuberculocidal

Packing : 100 ml, 250 ml , 500 ml , 1 lt, 5 lt



Ref.No	Product Name	Ship. Type
DC2060-100	Detro Derm Lotion	100ml
DC2060-500	Detro Derm Lotion	500ml
DC2060-1	Detro Derm Lotion	1lt
DC2060-5	Detro Derm Lotion	5lt
DC2062	Detro Derm Lotion Stand	



Superox Antisepsis

Disinfection Solution

Superox contains 6,2 pH and 200 ppm active chlorine (Hypochlorous Acid), 99.98% water, 0.02% HOCl. The product is used in humans as hand-foot, skin hygiene, mouth, nose, cleansing and mouthwash. The product is hand-foot, skin hygiene, mouth, nose, cleansing and gargle solution, skin around the eyes and hygiene.

Medical, public and personal conditions; Surgical hand washing, rubbing disinfectant is used directly. Gargle solution is used for skin cleansing and hygiene around the eyes, as post-surgical dressing solutions, for cleaning and moisturizing the wound area, it is used directly by spraying, pouring or wiping the area to be applied.

Cells called white blood cells in our body are for a physiological substance called hypochlorous Acid (HOCl). This substance is the largest and most effective spectrum known in nature, biocide. Thanks to its weak acidic feature, it is not irritating and harmful to human health. Effective is the best antiviral and antibacterial disinfectant known in nature, and the human body produces the best known disinfectant.

Certificates : CE,ISO 45001

Advantages:

- Natural antiseptic effect
- Has been proven by scientific studies
- Skin friendly
- Has most suitable pH value

Microbiological Impact:Bactericidal (including TBC), virucidal, fungicidal

Content: Water, Hypochlorous acid

Packaging : 100 ml. - 250 ml. - 500 ml. - 1 lt. - 5 lt.



Ref.No	Product Name	Ship. Type
DC2080-100	Superox Antisepsis	100ml
DC2080-250	Superox Antisepsis	250ml
DC2080-500	Superox Antisepsis	500ml
DC2080-1	Superox Antisepsis	1lt
DC2080-5	Superox Antisepsis	5lt



Copyright © 2020 Öncü Dental A.Ş. | All Rights Reserved.

Our registered brands.

T-protect | UltraFace® | FoldaFace® | T-Mask® | BodySafe® | Provisor Mask® | GoggleFace® | ProFeet®

dentac[®]

www.dentac.com | info@dentac.com

